

UNIVERSIDAD AUTÓNOMA DE MADRID

COMPUTER SCIENCE DEPARTMENT

---

# Biomedical Image Processing

## Assignment - 2

---

Roberto MARABINI RUIZ

## Contents

<b>1</b>	<b>Goal</b>	<b>2</b>
<b>2</b>	<b>Introduction</b>	<b>2</b>

# 1 Goal

To better understand the contrast transfer function in electron microscopy.

## 2 Introduction

The contrast transfer function describes how an object examined in a transmission electron microscope is imaged, essentially providing a description of distortions due to imperfect image formation by the microscope. By considering the recorded image as a CTF-degraded true object, describing the CTF allows the true object to be reverse-engineered.

In the following we introduce a first order approximation used to describe the CTF.

$$\begin{aligned}CTF(\mathbf{R}) &= w \cos(\gamma(\mathbf{R})) - \sqrt{1 - w^2} \sin(\gamma(\mathbf{R})) \\ \gamma(\mathbf{R}) &= 180\lambda(-\Delta Z\|\mathbf{R}\|^2 + \frac{Cs10^6\|\mathbf{R}\|^4\lambda^2}{2}) \text{ in degrees}\end{aligned}$$

where  $\mathbf{R}$  is the spatial frequency,  $\Delta Z$  denotes defocus in nm.,  $w$  is the percentage of amplitude contrast and  $Cs$  is the spherical aberration in mm. Finally, the factor  $10^6$  converts  $Cs$  from mm to nm. and  $\gamma(\mathbf{R})$  is termed in the specialized literature as wave aberration function.

Questions 1:

- Using matlab plot the CTF using the following two sets of values:
  - $Cs=0.6$  mm. accelerating voltage 300keV ( $\lambda=1.97$  pm),  $\Delta Z = 39.7$  nm.
  - $Cs=0.6$  mm. accelerating voltage 300keV ( $\lambda=1.97$  pm),  $\Delta Z = 90$  nm.
- For which frequencies the contrast is inverted?

## Questions 2:

- Simulate the aberrations introduced by the microscope by convolving the Einstein image with the above CTF.
- restore partially the image by multiplying (in Fourier space) it by the sign of the CTF
- restore partially the image by multiplying (in Fourier space) it by the  $\frac{1}{CTF+\epsilon}$  where  $\epsilon$  is a small number that avoid division by zero

## Questions 3:

- Create a uniform noise image and convolve it with the CTF. Plot the magnitude of the Fourier transform, the rings you see are the so-called Thon rings. (Fourier transform the image, multiply it by the CTF, plot the Fourier transform magnitude)
- What happen if you increase  $\Delta Z$ ? Does the distance between rings increase or decrease with a defocus increment?